Plastron 1/09

Tél/Fax: + 33 3 89 79 54 27

@: jkeilertplastron@orange.fr

## **Data Sheet**

## PLASTSTATIC 25 PE

**Application:** Internal antistatic for polyolefin; it gives almost instantaneous antistatic performance.

However in a lot of cases the effect might be only temporarily.

Applicable for PE; PP

Active Ingredients: glycerol ester

Active Content: 25%

**Carrier:** Modified polyethylene

Physiology: EEC No limitations for all polymers

USA No limitations for all polymers

Supplied as: whitish granules

## Dosage:

 Polyolefin films:
 < 50 μm</td>
 0.6 – 1.0 %

 > 50 μm
 1.0 – 2.0 %

Polyolefin moulded articles: 2.0 – 4.0 %

**Storage:** It is recommended to store the PLASTSTATIC 25-PE preferably below 30℃

The information given in this data sheet is based on many years of practical experience. Nevertheless can it only be a recommendation and a basis for practical testing, where the optimum dosage has to be determined individually.

All information given herein is only a recommendation. Its use shall remain the sole responsibility of the customer. No patent guarantee whatsoever shall be implied.